


# BeAT-HF Clinical Trial Summary


- Pre-market approval in 2019 based on symptomatic improvements and reductions in NT-proBNP
- Sustained and significant symptomatic improvements and safety at long-term follow-up
- Reductions in all-cause death, LVAD or transplant
- Improvement in the hierarchical composite (win ratio) of CV mortality, HF morbidity and QOL
- No statistically significant difference in the primary endpoint of CV mortality and HF morbidity



## Trial design

Design	Eligibility criteria
<ul style="list-style-type: none"> <li>• Prospective, multicenter, randomized 2-arm parallel-group, open-label with blinded evaluation</li> <li>• 73 US centers and 1 United Kingdom center</li> <li>• Randomized arms                             <ul style="list-style-type: none"> <li>• Barostim (BAT) plus GDMT (Barostim arm)</li> <li>• GDMT alone (Medical management arm)</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>• NYHA Class III or Class II (with a recent history of Class III)</li> <li>• Left ventricular EF <math>\leq</math> 35%</li> <li>• 6MHW 150 – 400 m</li> <li>• HF Hospitalization or NT-proBNP &gt; 400</li> <li>• Stable optimal management <math>\geq</math> 4 weeks</li> <li>• No class I indication for CRT</li> <li>• NT-proBNP &lt; 1600 pg/ml</li> </ul>

## BeAT-HF two-phase clinical trial



# Baseline demographics

	Barostim (n=163)	Control (n=160)
<b>Demographics</b>		
Age at Screening (years)	63 ± 11	63 ± 10
Gender (Female)	28 (17.2%)	35 (21.9%)
Race (Caucasian)	120 (73.6%)	116 (72.5%)
<b>Heart failure and physical status</b>		
SBP (mmHg)	120 ± 16	121 ± 16
DBP (mmHg)	74 ± 10	73 ± 10
HR (bpm)	75 ± 10	75 ± 11
BMI (kg/m <sup>2</sup> )	31 ± 5	31 ± 5
eGFR	62.5 ± 16.3	61.1 ± 18.9
NYHA: Class III	155 (95.1%)	151 (94.4%)
LVEF (%)	27 ± 6	28 ± 6
6 Minute Walk (m)	314 ± 66	300 ± 71
QOL	53 ± 24	51 ± 24
NT-proBNP (pg/mL)	736 (474, 1057)	704 (442, 1044)
LBBB	4 (2.5%)	2 (1.3%)
>=1 HF Hospitalization	66 (40.5%)	79 (49.4%)
Number of HF Hospitalizations	0.6 ± 0.9	0.7 ± 0.8

	Barostim (n=163)	Control (n=160)
<b>Co-Morbidities</b>		
Coronary Artery Disease	104 (63.8%)	107 (66.9%)
Atrial Fibrillation	53 (32.5%)	66 (41.3%)
Stroke or TIA	29 (17.8%)	37 (23.1%)
Chronic Kidney Disease	45 (27.6%)	43 (26.9%)
Type II Diabetes	74 (45.4%)	80 (50.0%)
<b>Heart failure treatment</b>		
Number of Meds	4.0 ± 1.3	4.1 ± 1.5
ACE-I / ARB / ARNI	143 (88%)	129 (81%)
ARNI	57 (35%)	43 (27%)
Beta-Blocker	152 (93%)	147 (92%)
MRA	74 (45%)	64 (40%)
Diuretic	138 (85%)	139 (87%)
Ivabradine	4 (2.5%)	9 (5.6%)
ICD	125 (77%)	127 (79%)


No significant difference between BAT and Control

## Sustained symptomatic improvement & safety


### Exercise capacity (6MHW)<sup>2</sup>



### Quality of life (MLWHF)<sup>2</sup>



### Functional status (NYHA Class)<sup>2</sup>




### Safety profile: MANCE<sup>2</sup>

97%

Freedom from major adverse neurological or cardiovascular system or procedure-related event rate in the Barostim arm


# Primary composite endpoint: CV mortality & HF morbidity<sup>5</sup>



No statistically significant difference between BAT and Control

\*Cumulative events per patient and rate ratio (RR) of treatment / control and 95% confidence interval estimated by negative binomial method

## Freedom from all-cause death, LVAD, and transplant<sup>2</sup>



-34%  
Relative reduction


Patients in the Barostim arm had a 34% reduction in all-cause death or the use of LVAD or heart transplant at > 4 years follow-up

\*Curves estimated using Kaplan-Meier method. Hazard ratio and p-value from Cox proportional hazards model.

## Hierarchical composite using win ratio analysis<sup>2</sup>

### Rationale:

- CV Mortality + HF Morbidity: 40% of patients contributed to the endpoint
- Win ratio: 100% of patients contribute to the endpoint




$$\text{Win Ratio} = \frac{\text{Total wins for treatment arm}}{\text{Total wins for control arm}}$$

1.26

HR 1.26 (95% CI 1.02, 1.58)

26% more patients in the Barostim arm had better outcomes as evaluated by hierarchical composite of CV mortality, LVAD/transplant, HF hospitalization and quality of life

# BeAT-HF summary of key evidence<sup>5</sup>



1. Zile MR, et al. J Am Coll Cardiol. 2020;76(1):1-13; 2. Instructions for Use 900133-001 Rev. D available at [www.cvr.com/ifu/](http://www.cvr.com/ifu/); 3. Gremeaux V, et al. Arch Phys Med Rehabil. 2011;92(4):611-619; 4. Rector TS, et al. J Card Fail. 1995;1(3):201-216; 5. Zile M, Presented at THT 2023, March 21, 2023.

## Important Safety Information

**CAUTION:** Federal law restricts this device to sale by or on the order of a physician. See System Reference Guide 900120-001 or 900133-001 for a complete instruction for use and a description of indications, contraindications, warnings, precautions and adverse events.

## Barostim™ Brief Summary for Physicians

The Barostim™ System is indicated for the improvement of symptoms of heart failure—quality of life, six-minute hall walk and functional status—for patients who remain symptomatic despite treatment with guideline-directed medical therapy, are NYHA Class III or Class II (who had a recent history of Class III), have a left ventricular ejection fraction  $\leq$  35%, a NT-proBNP < 1600 pg/ml and excluding patients indicated for Cardiac Resynchronization Therapy (CRT) according to AHA/ACC/ESC guidelines.

Patients are contraindicated if they have been assessed to have bilateral carotid bifurcations located above the level of the mandible, baroreflex failure or autonomic neuropathy, uncontrolled symptomatic cardiac bradyarrhythmias, carotid atherosclerosis that is determined by ultrasound or angiographic evaluation greater than 50%, ulcerative plaques in the carotid artery as determined by ultrasound or angiographic evaluation, known allergy to silicone or titanium.

**Warnings include:** only trained physicians may use this system, prescribing physicians should be experienced in the diagnosis and treatment of heart failure and should be familiar with the use of this system, monitor blood pressure and heart rate during Carotid Sinus Lead placement and when adjusting stimulation parameters intra-operatively, post-implantation, program the system to avoid the following: heart rate falls below 50 beats per minute (BPM), or systolic pressure falls below 90 mmHg, or diastolic blood pressure falls below 50 mmHg, or problematic adjacent tissue stimulation is noted, or undesirable interaction indicated by monitoring of any other implanted electrical device (see "Device Interaction Testing" in Section 10), or any other potentially hazardous patient responses are observed. Improper system implantation could result in serious injury or death. Do not use diathermy therapy including shortwave, microwave, or therapeutic ultrasound diathermy on patients implanted with the system. Patients should be counseled to stay at least 15 cm (6 inches) away from devices with strong electrical or magnetic fields such as strong magnets, loudspeaker magnets, Electronic Article Surveillance (EAS) system tag deactivators, arc welders, induction furnaces, and other similar electrical or electromechanical devices. This would include not placing items such as earphones in close proximity to the implanted pulse generator. The IPG may affect the operation of other implanted devices such as cardiac defibrillators, pacemakers, or neurological stimulation systems. For patients who currently have an implanted electrical medical device, physicians must verify compatibility with the implanted device during implantation of the system. Contralateral implant of the Barostim System may help to reduce potential interactions. Interactions are more likely in devices that contain a sensing function, such as an implantable cardiac defibrillator or pacemaker. If an interaction is observed, the Barostim System should be programmed to reduced therapy output settings in order to eliminate the interaction. If necessary, change settings in the other implant only if the changes are not expected to negatively impact its ability to perform its prescribed therapy. During the implant procedure, if device interactions cannot be eliminated the Barostim System should not be implanted.

**Precautions include:** the system should be implanted and programmed carefully to avoid stimulation of tissues near the electrode or in the area of the IPG pocket. Such extraneous stimulation could involve the following: the regional nerves, causing laryngeal irritation, difficulty swallowing, or dyspnea, the cervical musculature, causing intermittent contraction, skeletal muscles, causing intermittent contraction around the IPG pocket. Proper sterile technique during implantation should be practiced and aggressive pre-operative antibiotics are recommended. Infections related to any implanted device are difficult to treat and may necessitate device explantation.

It is anticipated that subjects will be exposed to operative and post-operative risks similar to related surgical procedures involving the neck and/or a pacemaker implant. These risks and potential risks of chronic device based baroreflex activation may include, but are not limited to: stroke, transient ischemic attack (TIA), systemic embolization, surgical or anesthetic complications, infection, wound complications, arterial damage, pain, nerve damage/stimulation, hypotension, hypertensive crisis, respiratory, exacerbation of heart failure, cardiac arrhythmias, tissue erosion/IPG migration, injury to baroreceptors, fibrosis, allergic reaction, general injury to user or patient, need for reoperation, secondary operative procedure, and death. Patients implanted with the system may receive Magnetic Resonance Imaging (MRI) only when all MR Conditional safety parameters are met as listed in the instructions for use.

CVRx, Barostim, BAT, Barostim NEO, Barostim NEO2, BATwire, Outsmart the heart, and Barohub are all trademarks of CVRx, Inc. © 2024 CVRx, Inc. All rights reserved. All other trademarks are property of their respective owners.

For a list of all potential benefits and risks go to [www.cvr.com/benefit-risk-analysis/](http://www.cvr.com/benefit-risk-analysis/)

CVRx, Inc.  
9201 West Broadway Avenue, Suite 650  
Minneapolis, MN 55445 USA  
P 763-416-2840 | F 763-416-2841  
[cvrx.com](http://cvrx.com)

**CVRx**  
Outsmart the heart